CareFusion

Smoke Check Meter

Operating Manual – English

Operating Manual – CareFusion USA

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Contents

I.	Overview.....	3
	Introduction – English.....	4
	Cigarette Consumption	4
	Package Contents.....	5
	Operation.....	7
	Calibration	10
	Battery Life	12
	Sensor Expiry	13
	Cleaning	13
	Servicing.....	13
	Symbols.....	14
	Environment.....	14
	Specifications	14
II.	Overview.....	17
	Introduction – CareFusion USA	17
	Cigarette Consumption	18
	Package Contents.....	19
	Operation.....	21
	Calibration	24
	Battery Life	26
	Sensor Expiry	27
	Cleaning	27
	Servicing.....	27
	Symbols.....	28
	Environment.....	28
	Specifications	28
III.	Présentation	30
	Introduction – Français	31
	Consommation de Cigarettes	31

Contenu de l’emballage.....	32
Instructions	34
Calibration	36
Durée de vie pile.....	40
Expiration de la Cellule	40
Nettoyage	40
Informations techniques.....	40
Symbols.....	40
Environment.....	41
Spécifications	41
 IV. Überblick.....	 42
Einführung – Deutsch	43
Zigarettenkonsum	44
Packungsinhalte	45
Bedienung.....	47
Eichung.....	50
Lebensdauer der Batterie	53
Verfall des Sensors.....	54
Reinigung.....	54
Wartung.....	55
Symbole	55
Umgebung	55
Anforderungen.....	55
 Customer contact information	 57

Smoke Check Meter

Operating manual

Overview

The Smoke Check is a hand held battery operated device used to measure the concentration of carbon monoxide, CO, on the breath.

It is used for instructional purposes in smoking cessation programs and to check on the clients' progress and compliance.

It is accurate, easy to use, and has many features designed to simplify its operation.

These include:

- Auto zero function
- Smoking level colour light indicators
- High CO alarm
- Simple calibration

Please note: It is recommended that this unit be calibrated upon receipt

I. Introduction - English

The Smoke Check is based on an electrochemical fuel cell sensor, which works through the reaction of carbon monoxide (CO) with an electrolyte at one electrode, and oxygen (from ambient air) at the other. This reaction generates an electrical current proportional to CO concentration. Output from the sensor is monitored by a microprocessor, which detects and displays peak expired concentrations of alveolar gas. High levels of expired CO indicate raised levels of carboxyhaemoglobin, most commonly caused by cigarette smoking.

The results are displayed in four ranges on a clear LCD display. Warning lights are provided to give an instant indication of the smoking level. A particular feature of the CareFusion Smoke Check is its stability of calibration and low cross sensitivity to other exhaled gases such as hydrogen.

References

Jarvis MJ, Belcher M, Vesey C, Hutchison DCS

Low cost carbon monoxide monitors in smoking assessment

Thorax 1986; 41: 886-887

Cigarette Consumption

The Smoke Check meter provides a simple screening test for cigarette consumption for use in anti-smoking clinics and all smoking cessation programs.

Measurement of expired CO has been well validated as an indirect measure of cigarette consumption and is widely used in smoking cessation programs.

Typical values for expired CO in smokers, together with the alarm light status, are given below:

CO(ppm)	Cigarette consumption	Indicator
0 - 6	Non smoker	Green
7 -10	Light smoker	Amber
11 – 20	Heavy smoker	Red
20+	Very heavy smoker	Flashing red + alarm

Please note that some urban areas may have high environmental levels of CO. This can cause a rise in exhaled CO of a few ppm above that which is normally present on the breath. In these cases it is possible for a non-smoker to appear at the bottom of the 'light smoker' range (7 – 10 ppm).

References

Jones RH, Ellicott MF, Cadigan JB, Gaensler EA

The relationship between alveolar and blood carbon monoxide concentrations during breath holding

Journal of Laboratory and Clinical Medicine 1958; 51, 553 - 564

Wald NJ, Idle M, Boreham J, Baily A





Carbon monoxide in breath in relation to smoking and carboxyhaemoglobin levels

Thorax 1981; 36, 366-369

Package Contents

The Smoke Check is supplied with a carrying case containing the following items:

1. Smoke Check meter (Cat No. 36-SC01-STK)
2. 22 mm mouthpiece adapter (Box of 10 Cat No.36-PSA2000)
3. 22 mm reducing connector (Cat No.36-MEC1007)
4. PP3 Battery (Cat No.36-BAT1002)
5. Disposable mouthpieces (Box 250 Cat No.36-PSA1200)
6. Calibration button tool (Cat No. 36-MEC1299)


Operation

For accurate results the Smoke Check should be used at room temperature. If the instrument has been stored in cool or hot conditions then allow time to reach room temperature prior to use.

Install the PP3 battery by sliding open the battery cover, clipping the battery in place and replacing the cover. Insert the mouthpiece adapter into the Smoke Check meter and then insert a disposable cardboard mouthpiece into the adapter. If possible, the subject should rinse their mouth with clean water prior to performing a test. Do not use mouthwash, as the readings will be affected.

To obtain an accurate reading from a sample of alveolar gas the patient should inspire and hold their breath for 20 seconds before expiring slowly and fully through the mouthpiece. As an aid to timing the breath holding period, a buzzer will sound 20 seconds after the unit is turned on. However the instrument may be used any time after the first 12 seconds when the auto-zero has been performed.

Turn the unit on by pushing the central slide switch up and the three coloured lights will illuminate momentarily together with all the display legends:

Bat	7-10 ppm	20+ ppm
Wait		
Blow	0-6 ppm	11-20 ppm
Gas		


Ask the subject to inspire and hold their breath.

The unit will auto-zero for twelve seconds during which time the following will be displayed:


Important note: the unit must not be used during this period.

After twelve seconds the display will change to:

The subject may now seal their lips around the mouthpiece and exhale slowly and fully if they cannot hold their breath any longer. However, more accurate results will be obtained if they hold their breath until the buzzer sounds and the display changes to:

The CO reading will rise to a plateau over the course of several seconds. The final value will be held until the unit is turned off and will be displayed in one of the four ranges:


- 0–6 ppm with green indicator
- 7–10 ppm with amber indicator
- 11–20 ppm with red indicator
- 20+ ppm with flashing red indicator and alarm

Important note:

Before repeating a measurement the unit must be turned off, and the mouthpiece and adapter removed for at least 1 minute. This is to allow re-equilibration with ambient air and to dry the surface of the sensor. Visually inspect that all moisture has evaporated from the surface of the sensor before reuse.

If the unit is switched on again too quickly after use there may be a response to residual expired carbon monoxide from the previous test.

In this case the display will show:

If this is displayed then turn the unit off, remove the mouthpiece adapter, and expose to ambient air for 2 minutes before repeating the test.

Note: If this warning appears again after following the above procedure then turn off and leave the sensor exposed to ambient air for a further 3 minutes. If the same message appears again then this indicates possible contamination of the sensor with a solvent.

In this case remove all sources of solvent from the vicinity of the sensor and expose to ambient air for 24 hours before switching on again.


Calibration

Calibration will remain stable to within 2% over one month and typically to within 10% over 6 months. CareFusion supplies calibration gas (20ppm CO in air) complete with a control valve for simple and economic user calibration.

The gas is supplied in convenient, disposable, aluminium aerosol cans containing 18 litres of gas (Cat. No.36-MCGK020).


Each cylinder is capable of completing 10 calibration procedures.

Locate the calibration button positioned inside a hole on the right hand side of the instrument as shown below.

Turn the unit on and wait for the unit to display 'Blow'.

Screw the control valve firmly onto the cylinder and connect the gas supply as shown below:

The plastic tubing supplied with the gas should be pushed firmly over the reducing connector.

Slowly turn the control knob anti-clockwise until the ball in the flow indicator is between the two marks. This will then supply a gas flow of approximately 0.25 l/min. Apply this flow for 25 seconds and then use the calibration button tool to depress the calibration button.

The sounder will beep and the display will flash three times as the calibration is recorded. The 'GAS' legend will then be displayed. Turn off the gas supply and the unit. When the cell has expired it will be impossible to obtain a correct calibration and all the legends and the three lights will flash five times with no beeps after the calibration button is depressed. The 'GAS' legend will be displayed as before and the unit must be turned off before recalibration is attempted.


Important Notes:

- Only certified calibration gas (20ppm) from a reputable source should be used.
- Ensure that no CO is present on the sensor for 3 minutes before starting the calibration procedure.
- Ensure that the instrument and gas cylinder have stabilised at room temperature before calibrating.
- Do not immerse the unit in the calibration gas.
- In order to stop any gas leak from the can after calibration, ensure that the control knob is tightened firmly.

Battery Life

Battery life is sufficient for approximately 8000 tests.

When the battery has approximately 1 hour of useful life left the following message will be displayed:

and the alarm will sound momentarily after the unit is first switched on. When the battery has completely expired the above will be displayed continuously and the battery must be replaced.

Note: Please remove the battery if the meter is likely to be unused for some time.

Sensor Expiry

Sensor life ranges from 2 to 5 years and depends upon both the amount of exposure to CO and other gases, particularly solvents such as alcohol and cleaning fluids.

When the sensor has expired it will become impossible to obtain a correct calibration. When this occurs, the Smoke Check must be returned to CareFusion or an authorised agent, for sensor replacement.

Cleaning

The mouthpiece adapter may be cleaned using a mild detergent solution or replaced if there is risk of contamination. Exposed surfaces of the meter, other than the sensor surface, may be wiped with a clean, damp cloth. It is recommended that this procedure is performed after each use and all used cardboard mouthpieces discarded.

Important note: The sensor surface must not be wiped with any aqueous solutions and must not be exposed to solvents e.g. alcohol or permanent damage may result.

Servicing

A full service manual including circuit diagram and parts list is available upon request.

Symbols

Type B device

0086

In accordance with Directive 93/42/EEC

Disposal in compliance with WEEE

Environment

This instrument complies with directive EN60601-1-2 electromagnetic compatibility but can be affected by cellular phones and by electromagnetic interference exceeding levels specified in EN 50082-1:1992

Specifications

Sensor type	Electro-chemical fuel cell
Range	0 - 20 ppm
Detected levels:	Display:
0 to 6ppm	0-6 ppm with green indicator
7 to 10ppm	7-10 ppm with amber indicator
11 to 20ppm	11-20 ppm with red indicator
>20 ppm	20+ ppm with red indicator + alarm
Accuracy	+/-5% of reading or 1ppm whichever is the greater
Sensitivity drift	0.5%/°C
Sensor life	2 to 5 years
Response time	< 15 sec (to 90% of reading)
Hydrogen cross sensitivity	<15%
Operating temperature	15 – 25 °C
Operating pressure	Atmospheric +/- 10%

Pressure coefficient	0.02% signal per mBar
Relative humidity (Non condensing)	15 - 90% continuous (0 - 99% intermittent)
Baseline drift	0ppm (auto-zero)
Long term drift	< 2% signal loss per month
Power source	Single Alkaline 9 volt PP3
Battery life	> 8000 tests
Weight	130g without battery
Dimensions	170 x 60 x 26 mm
Display	Custom LCD
Storage temperature	-20 ° to +70 °C
Storage humidity	30% to 90%

Smoke Check Meter

Operating manual

Federal (USA) law restricts this device to sale by or on the order of a physician or licensed practitioner.

CareFusion,
22745 Savi Ranch Parkway,
Yorba Linda
CA 92887-4668
USA

Overview

The Smoke Check is a hand held battery operated device used to measure the concentration of carbon monoxide, CO, on the breath.

It is used for instructional purposes in smoking cessation programs and to check on the clients' progress and compliance.

It is accurate, easy to use, and has many features designed to simplify its operation.

These include: Auto zero function
Smoking level colour light indicators
High CO alarm
Simple calibration

Please note: It is recommended that this unit be calibrated upon receipt

II. Introduction – CareFusion USA

The Smoke Check is based on an electrochemical fuel cell sensor, which works through the reaction of carbon monoxide (CO) with an electrolyte at one electrode, and oxygen (from ambient air) at the other. This reaction generates an electrical current proportional to CO concentration. Output from the sensor is monitored by a microprocessor, which detects and displays peak expired concentrations of alveolar gas. High levels of expired CO indicate raised levels of carboxyhaemoglobin, most commonly caused by cigarette smoking.

The results are displayed in four ranges on a clear LCD display. Warning lights are provided to give an instant indication of the smoking level. A particular feature of the Smoke Check is its stability

of calibration and low cross sensitivity to other exhaled gases such as hydrogen.

References

Jarvis MJ, Belcher M, Vesey C, Hutchison DCS

Low cost carbon monoxide monitors in smoking assessment

Thorax 1986; 41: 886-887

Cigarette Consumption

The Smoke Check meter provides a simple screening test for cigarette consumption for use in anti-smoking clinics and all smoking cessation programs.

Measurement of expired CO has been well validated as an indirect measure of cigarette consumption and is widely used in smoking cessation programs.

Typical values for expired CO in smokers, together with the alarm light status, are given below:

CO(ppm)	Cigarette consumption	Indicator
0 - 6	Non smoker	Green
7 -10	Light smoker	Amber
11 – 20	Heavy smoker	Red
20+	Very heavy smoker	Flashing red + alarm

Please note that some urban areas may have high environmental levels of CO. This can cause a rise in exhaled CO of a few ppm above that which is normally present on the breath. In these cases it is possible for a non-smoker to appear at the bottom of the 'light smoker' range (7 – 10 ppm).

References

Jones RH, Ellicott MF, Cadigan JB, Gaensler EA

The relationship between alveolar and blood carbon monoxide concentrations during breath holding

Journal of Laboratory and Clinical Medicine 1958; 51, 553 - 564

Wald NJ, Idle M, Boreham J, Baily A





Carbon monoxide in breath in relation to smoking and carboxyhaemoglobin levels

Thorax 1981; 36, 366-369

Package Contents

The Smoke Check is supplied with a carrying case containing the following items:

1. Smoke Check meter (Cat No. 36-SC01-STK)
2. 22 mm mouthpiece adapter (Box of 10 Cat No.36-PSA2000)
3. 22 mm reducing connector for calibration (Cat No.36-MEC1007)
4. PP3 Battery (Cat No.36-BAT1002)
5. Disposable mouthpieces (Box 250 Cat No.36-PSA1200)
6. Calibration button tool (Cat No. 36-MEC1299)

Operation

For accurate results the Smoke Check should be used at room temperature. If the instrument has been stored in cool or hot conditions then allow time to reach room temperature prior to use.

Install the PP3 battery by sliding open the battery cover, clipping the battery in place and replacing the cover. Insert the mouthpiece adapter into the Smoke Check meter and then insert a disposable cardboard mouthpiece into the adapter. If possible, the subject should rinse their mouth with clean water prior to performing a test. Do not use mouthwash, as the readings will be affected.


To obtain an accurate reading from a sample of alveolar gas the patient should inspire and hold their breath for 20 seconds before expiring slowly and fully through the mouthpiece. As an aid to timing the breath holding period, a buzzer will sound 20 seconds after the unit is turned on. However the instrument may be used any time after the first 12 seconds when the auto-zero has been performed.

Turn the unit on by pushing the central slide switch up and the three coloured lights will illuminate momentarily together with all the display legends:

Bat	7-10 ppm	20+ ppm
Wait		
Blow	0-6 ppm	11-20 ppm
Gas		


Ask the subject to inspire and hold their breath.

The unit will auto-zero for twelve seconds during which time the following will be displayed:


Important note: the unit must not be used during this period.

After twelve seconds the display will change to:

The subject may now seal their lips around the mouthpiece and exhale slowly and fully if they cannot hold their breath any longer. However, more accurate results will be obtained if they hold their breath until the buzzer sounds and the display changes to:

The CO reading will rise to a plateau over the course of several seconds. The final value will be held until the unit is turned off and will be displayed in one of the four ranges:


- 0–6 ppm with green indicator
- 7–10 ppm with amber indicator
- 11–20 ppm with red indicator
- 20+ ppm with flashing red indicator and alarm

Important note:

Before repeating a measurement the unit must be turned off, and the mouthpiece and adapter removed for at least 1 minute. This is to allow re-equilibration with ambient air and to dry the surface of the sensor. Visually inspect that all moisture has evaporated from the surface of the sensor before reuse.

If the unit is switched on again too quickly after use there may be a response to residual expired carbon monoxide from the previous test.

In this case the display will show:

If this is displayed then turn the unit off, remove the mouthpiece adapter, and expose to ambient air for 2 minutes before repeating the test.

Note: If this warning appears again after following the above procedure then turn off and leave the sensor exposed to ambient air for a further 3 minutes. If the same message appears again then this indicates possible contamination of the sensor with a solvent.

In this case remove all sources of solvent from the vicinity of the sensor and expose to ambient air for 24 hours before switching on again.


Calibration

Calibration will remain stable to within 2% over one month and typically to within 10% over 6 months. CareFusion supplies calibration gas (20ppm CO in air) complete with a control valve for simple and economic user calibration.

The gas is supplied in convenient, disposable, aluminium aerosol cans containing 18 litres of gas.


Each cylinder is capable of completing 10 calibration procedures.

Locate the calibration button positioned inside a hole on the right hand side of the instrument as shown below.

Turn the unit on and wait for the unit to display 'Blow'.

Screw the control valve firmly onto the cylinder and connect the gas supply as shown below:

The plastic tubing supplied with the gas should be pushed firmly over the reducing connector.

Slowly turn the control knob anti-clockwise until the ball in the flow indicator is between the two marks. This will then supply a gas flow of approximately 0.25 l/min. Apply this flow for 25 seconds and then use the calibration button tool to depress the calibration button.

The sounder will beep and the display will flash three times as the calibration is recorded. The 'GAS' legend will then be displayed. Turn off the gas supply and the unit. When the cell has expired it will be impossible to obtain a correct calibration and all the legends and the three lights will flash five times with no beeps after the calibration button is depressed. The 'GAS' legend will be displayed as before and the unit must be turned off before recalibration is attempted.


Important Notes:

- Only certified calibration gas (20ppm) from a reputable source should be used.
- Ensure that no CO is present on the sensor for 3 minutes before starting the calibration procedure.
- Ensure that the instrument and gas cylinder have stabilised at room temperature before calibrating.
- Do not immerse the unit in the calibration gas.
- In order to stop any gas leak from the can after calibration, ensure that the control knob is tightened firmly.

Battery Life

Battery life is sufficient for approximately 8000 tests.

When the battery has approximately 1 hour of useful life left the following message will be displayed:

and the alarm will sound momentarily after the unit is first switched on. When the battery has completely expired the above will be displayed continuously and the battery must be replaced.

Note: Please remove the battery if the meter is likely to be unused for some time.

Sensor Expiry

Sensor life ranges from 2 to 5 years and depends upon both the amount of exposure to CO and other gases, particularly solvents such as alcohol and cleaning fluids.

When the sensor has expired it will become impossible to obtain a correct calibration. When this occurs, the Smoke Check must be returned to CareFusion, or an authorised agent, for sensor replacement.

Cleaning

The mouthpiece adapter may be cleaned using a mild detergent solution or replaced if there is risk of contamination. Exposed surfaces of the meter, other than the sensor surface, may be wiped with a clean, damp cloth. It is recommended that this procedure is performed after each use and all used cardboard mouthpieces discarded.

Important note: The sensor surface must not be wiped with any aqueous solutions and must not be exposed to solvents e.g. alcohol or permanent damage may result.

Servicing

A full service manual including circuit diagram and parts list is available upon request.

For US Customer Care: Toll Free 1-800-231-2466, phone 1-714-283-2228

Symbols

Type B device

Environment

This instrument complies with directive EN60601-1-2 electromagnetic compatibility but can be affected by cellular phones and by electromagnetic interference exceeding levels specified in EN 50082-1:1992

Specifications

Sensor type	Electro-chemical fuel cell
Range	0 - 20 ppm
Detected levels:	Display:
0 to 6ppm	0-6 ppm with green indicator
7 to 10ppm	7-10 ppm with amber indicator
11 to 20ppm	11-20 ppm with red indicator
>20 ppm	20+ ppm with red indicator + alarm
Accuracy	+/-5% of reading or 1ppm whichever is the greater
Sensitivity drift	0.5%/°C
Sensor life	2 to 5 years
Response time	< 15 sec (to 90% of reading)
Hydrogen cross sensitivity	<15%
Operating temperature	15 – 25 °C
Operating pressure	Atmospheric +/- 10%
Pressure coefficient	0.02% signal per mBar
Relative humidity (Non condensing)	15 - 90% continuous (0 - 99% intermittent)
Baseline drift	0ppm (auto-zero)
Long term drift	< 2% signal loss per month

Power source	Single Alkaline 9 volt PP3
Battery life	> 8000 tests
Weight	130g without battery
Dimensions	170 x 60 x 26 mm
Display	Custom LCD
Storage temperature	-20 ° to +70 °C
Storage humidity	30% to 90%

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Customer contact information

UK Customers only

For all Sales Order processing for products, training and spare parts, Service and Technical Support enquiries, please contact the following:

CareFusion UK 232 Ltd
UK Customer Service & Support
The Crescent
Jays Close
Basingstoke
RG22 4BS

Customer Service Sales Enquiries:

Telephone: 01256 388550

Email: micro.uksales@carefusion.com

Factory Repair and Administration Enquiries:

Telephone: 01256 388552

Email: micro.ukservice@carefusion.com

Technical Support Enquiries:

Telephone: 01256 388551

Email: support.rt.eu@carefusion.com

International customers only

For all Sales Order processing for products and Spare parts, Service and Technical Support enquiries, please contact the following:

Carefusion Germany 234 GmbH
Customer Service & Support International
Leibnizstrasse 7
97204 Hoechberg
Germany

Customer Service Sales Enquiries:

Telephone: 0049 931 4972 670

Email: micro.internationalsales@carefusion.com

Factory Repair and Administration Enquiries:

Telephone: 0049 931 4972 867

Email: support.admin.eu@carefusion.com

Technical Support Enquiries:

Telephone: 0049 931 4972 127

Email: support.rt.eu@carefusion.com

Notes:

CareFusion UK 232 Ltd.,
Quayside
Chatham Maritime
Kent ME4 4QY
U.K.

JDE no. 36-MAN1281
Drg no. 041-04WW
Issue 1.0
January 2010
© CareFusion 2010